

### Varian 500 Tutorial for the acquisition of 1D $^1\text{H}$ and $^{13}\text{C}$ spectra

- 1) log in to the linux account.
- 2) double click on the VNMRJ icon to start the program.
- 3) check that the probe file is "OneNMR\_W036".
- 4) insert the sample tube in the spinner. Use the gauge to measure the depth.
- 5) click on < eject > in the "start/standard" tab, and take the sample out of the magnet.
- 6) put your sample on top of the magnet. Be sure that the spinner is floating before dropping it.
- 7) click <insert> in the "start/standard" tab.
- 8) the line commands "e" and "i" can also be used to eject and insert samples respectively.
- 9) wait until the sample is loaded in the magnet. The yellow light, in the switch box on top of the magnet, turns on when a sample is in the probe.
- 10) **Tune the probe** to  $^1\text{H}$ : type "protune" on the command line. The tuning window pops up.
- 11) click on < $^1\text{H}$ > and wait until the tuning is done. Then, close the tuning window.
- 12) **Lock to the solvent:** select the solvent on the tab "Start/standard".
- 13) on the "Start/Standard" tab, click on <Find z0> to lock the sample.

- 14) check whether there is a “Mapname:”, then click on <Gradient Autoshim>.
- 15) while doing gradient shims and Find Z0, the solvent will temporarily switch to D2O and spinning will turn off.
- 16) on the tab Start/ Spin/Temp, set the spinning rate to 20 Hz and then press the button <Regulate Speed>.
- 17) temperature can be set in the same tab. Select the temperature and then click on <Regulate Temp.>.
- 18) choose a 1D-<sup>1</sup>H experiment, click on the menu “experiment” and then select “proton”.
- 19) on the tab “Acquire” you can modify the experimental parameters as needed.
- 20) run the experiment, on the command line type “ga” and then press <enter>. Or just click on <acquire> in the “acquire” tab.
- 21) to display the spectrum, type “wft” and then press <enter>.
- 22) usefull commands: **wft**: to refresh the spectrum, **aph**: automatic phasing, **nt: number of scans vsadj** : adjust the high of the spectrum, **f**: brings back the full spectrum after it was zoomed. **movesw**: adjusts sw automatically, **movetof**: changes the central frequency.
- 23) saving files: in the file menu, select “save as”.

### **Printing spectra**

- 1) select the “process/plot” tab. Select the desire options and then click on <Auto Plot>.

### **Running a carbon experiment.**

- 1) **Tune <sup>13</sup>C**: type “protune” on the command line and then on the window that pops up, click on <<sup>13</sup>C>.

- 2) Remember that for a carbon experiment, you also need to tune  $^1\text{H}$ .
- 3) change to experiment 2 by typing “jexp2” and then press <enter>.
- 4) select the experiment “carbon”.
- 5) type on the command line “ga”, then press <enter>. Or just click on <acquire>.

**Once you are done, put the standard sample back in the magnet and lock it:**

- 1) select the solvent ( $\text{D}_2\text{O}$ ), then press <Find Z0>.
- 2) stop the spinning. On tab “start/ spin/temp” click on <spin off>.
- 3) set the temperature back to 26 °C.
- 4) close VNMRJ; select “File” on the upper left corner and then “Exit Vnmrj”.
- 5) log out from your account.
- 6) log in the log book.

### **Retreiving data from the spectrometer computer using SSH**

- 0) Open the WinSCP program.
  - 1) click on <NEW> .
  - 2) Fill in the windows that pops up: host name: juliet.chem.utk.edu, user name: .... , and password:.....
  - 3) click on <save>.
  - 4) click <OK> on next window.

- 5) when trying to connect for the first time, a windows pops up. Just click on <yes>.
- 6) the right panel is the host computer, change the directory to vnmrsys/data.
- 7) download the data.
- 8) once you are done, diconnect from the host computer.